

Operating Theatre Essentials

Incorporating

*Safe Sharps Management
& Tubing Essentials*

purple surgical

A message from the CEO...

Change is an inevitable part of human development and it's often the process of change that makes life most exciting.

Cory Bros (Surgical Instrument Makers) Ltd. was founded in England, in September 1909. Over the last 100 years, the Cory Bros name has become synonymous with high quality, cost effective surgical instrumentation and devices.

In 2000, Cory Bros created a new company, Purple Surgical International Ltd. to sell our products worldwide, where we are now represented in over 80 countries.

In 2002 our history turned full circle when we opened our current manufacturing facility in Taunton, Somerset and we are extremely proud to be one of the few remaining Great British surgical product manufacturers.

In 2012 we made the decision to coordinate our worldwide activities more closely, so from 2013, we shall cease using the Cory Bros name, with all our products now carrying the Purple Surgical brand.

In the 104 years since Cory Bros was founded, the company has become a highly respected part of the British healthcare market and we are looking forward to Purple Surgical becoming equally well known and loved.

As part of this re-branding exercise, we have chosen to combine our extensive range of Operating Theatre Essentials, Safe Sharps Management, and Tubing Essentials products into one brochure for the first time. This will enable you to more easily select the perfect combination of surgical accessories from one company, that best suit your requirements.

I very much hope that Purple Surgical can become a significant and valued partner in your day-to-day professional life.

Robert Sharpe
Chairman & CEO



Sharps Safety Station™

The unique, patented Purple Surgical Sharps Safety Station™ range provides the ultimate in sterile sharps organisation, counting and disposal of contaminated sutures, blades and needles. This extensive range incorporates a number of unique safety features, which significantly reduce the risk of inadvertent sharps injuries.

Features and benefits include:

- Patented, quick, simple and safe universal blade remover
- Unique, convenient integral scalpel shield for safe and simple scalpel storage and retrieval
- Unique locking mechanism prevents inadvertent opening whilst enabling re-counts if required
- Unique, secure, three hinge design prevents box opening during transportation
- Deep overlapping edges prevent accidental sharps migration
- Four size options to meet surgical requirements
- Wide choice of formats, including foam blocks, magnetic surfaces, adhesive surfaces and foam strips, to meet customer preference and clinical indication
- International bio-hazardous waste colour identification




Product	Code	Description	Size	Pack Qty
Sharps Safety Station™	PS2535	Double Magnet	Medium	25 / Box
Sharps Safety Station™	PS2560	Double Magnet	Large	25 / Box
Sharps Safety Station™	PS2635	Magnet and Foam Strips	Medium	25 / Box
Sharps Safety Station™	PS2615	Magnet and Foam Strips	Large	25 / Box
Sharps Safety Station™	PS2925	Magnet and Foam Block	Medium	25 / Box
Sharps Safety Station™	PS2940	Magnet and Foam Block	Large	25 / Box
Sharps Safety Station™	PS2630	Double Foam Strips	Medium	25 / Box
Sharps Safety Station™	PS2922	Adhesive and Foam Block	Medium	25 / Box

Sharps Count

The Purple Surgical Sharps Count is a rigid plastic, puncture-proof box for the organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:

- Strong rigid construction prevents accidental penetration by contaminated sharps
- A clear lid option allows convenient sharps visualisation and counting
- Unique, locking mechanism prevents inadvertent opening whilst enabling re-counts if required
- Unique, secure, three hinge design prevents box opening during transportation
- Deep overlapping edges prevent accidental sharps migration
- International bio-hazardous waste colour identification



Product	Code	Description	Size	Pack Qty
Sharps Count	PS3100	Magnet with Clear Top	Medium	25 / Box
Sharps Count	PS3110	Magnet with Clear Top	Large	25 / Box
Sharps Count	PS2200	Double Magnet	Medium	25 / Box
Sharps Count	PS2910	Magnet and Foam Block	Small	25 / Box

Puncture-Proof-Pad

The Purple Surgical Puncture-Proof-Pad is a unique rigid plastic, puncture-proof-pad for the simple, safe organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:

- Strong rigid construction prevents accidental penetration by contaminated sharps
- Wide adhesive strip ensures secure sharps retention
- Secure locking mechanism prevents inadvertent opening
- Overlapping edges prevent sharps migration
- An optional 'fold flat' edge allows the safe accommodation of hypodermic needles and scalpel handles when required
- International bio-hazardous waste colour identification



Product	Code	Description	Pack Qty
Puncture-Proof-Pad	PS2002	Adhesive and Foam	50 / Box

Purple Turtle

The Purple Turtle is a unique, cost effective puncture resistant pad for the simple, safe organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:

- Tough puncture resistant plastic offers excellent protection against the risk of accidental sharps penetration
- Wide adhesive strip ensures secure sharps retention
- Overlapping edges prevent sharps migration
- A castellated edge allows the safe accommodation of hypodermic needles



Product	Code	Description	Pack Qty
Purple Turtle	PS2001	Adhesive and Foam	60 / Box

Sharps Pad

The Purple Surgical Sharps Pad is a simple, low cost cardboard backed foam pad which offers the minimum requirement for the organisation, counting and disposal of contaminated sutures, blades and needles.

Features and benefits include:

- Purpose designed
- Convenient
- Low cost




Product	Code	Description	Pack Qty
Sharps Pad	PS2000	Adhesive and Foam	50 / Box

Scalpel Shield

The Purple Surgical Scalpel Shield is a purpose designed device for the safe, secure storage of up to three scalpels in the sterile field.

Features and benefits include:

- Conveniently stores handles in the 'upright' position ensuring safe and simple placement and retrieval
- Shrouded blade design prevents accidental personal injury and damage to drapes
- Adhesive tab secures and stabilises device in the sterile field



Product	Code	Description	Pack Qty
Scalpel Shield	PS1593	Scalpel Shield	25 / Box

Sharps Passing Tray

The Purple Surgical Sharps Passing Tray is a lightweight, purpose designed, single-use device that provides a secure 'hands free' transfer medium for passing sharp instruments, thus reducing the risk of third party inflicted sharps injuries.

Features and benefits include:

- Central cavity design ensures all sharp instruments are safely shielded during transfer from person to person
- Versatile cavity accepts scalpels, needle holders, syringes and all sharp instruments (e.g. scissors, trocars etc.)
- Two trays may be joined to allow safe, 'hands free' transfer of multiple instruments
- International bio-hazardous waste colour identification



Product	Code	Description	Pack Qty
Sharps Passing Tray	PS4594	Sharps Passing Tray	25 / Box

Mag-Mat

The Purple Surgical Mag-Mat is a reusable magnetic mat that provides a secure 'hands free' transfer zone for metallic instruments, thus reducing the risk of third party inflicted sharps injuries.

Features and benefits include:

- Strong ceramic magnetic surface ensures secure retention of metallic instruments
- Will not magnetise instruments
- Autoclavable
- High quality, reinforced mesh design ensures a cost effective, long lasting device
- Large (30cm x 40cm) surface allows for containment of multiple instruments
- Waterproof construction minimises the risk of any pathogen cross infection



Product	Code	Description	Pack Qty
Mag-Mat	PS4116	Reusable Mag-Mat 30cm x 40cm (12" x 16")	1 / Each

Magnetic Drapes

Purple Surgical Magnetic Drapes are lightweight, single-use and provide a secure 'hands free' transfer zone for metallic instruments, thus reducing the risk of third party inflicted sharps injuries.

Features and benefits include:

- Strong magnetic surface ensures secure retention of metallic instruments
- Convenient single-use device
- Lightweight construction makes the drape particularly suitable for young, elderly or frail patients
- Waterproof construction minimises the risk of any pathogen cross infection
- Available in two sizes



Product	Code	Description	Pack Qty
Magnetic Drape	PS4121	Single-use Magnetic Drape 25cm x 40cm (10" x 16")	25 / Box
Magnetic Drape	PS4122	Single-use Magnetic Drape 40cm x 50cm (16" x 20")	25 / Box

Golden Retriever

The Purple Surgical Golden Retriever is a reusable device for the quick, simple and safe retrieval of metallic sharps that have been inadvertently dropped onto the operating room floor.

Features and benefits include:

- Strong magnetic roller ensures safe, 'hands free' sharps retrieval
- Robust design ensures a cost effective, long lasting device
- Long handle enables retrieval without the need for bending



Product	Code	Description	Pack Qty
Golden Retriever	PS4000	Reusable Magnetic Retriever	1 / Each

Instrument Pouch

The Purple Surgical Instrument Pouch enables the secure, convenient housing of instruments close to the surgical site, thus minimising the necessity for hand-to-hand transfer.

Features and benefits include:

- Adhesive strip securely attaches the pouch to the patient drape
- Deep pockets ensure the safe retention of all standard length minimally invasive instruments
- Instrument separation is assured by the three pocket design
- The clear material design ensures quick instrument identification and retrieval



Product	Code	Description	Pack Qty
Instrument Pouch	PS4050	Instrument Pouch	25 / Box

Patient Eye Shields

The Purple Surgical Patient Eye Shields effectively cover and protect eyes from dust and dirt.

Features and benefits include:

- Universal design fits both left and right eyes, thus minimising inventory
- Transparent material allows limited vision if required
- High quality, smooth, ergonomic design for optimum patient comfort



Product	Code	Description	Pack Qty
Patient Eye Shield	PS5666	Universal Patient Eye Shield	100 / Box

Skin Marker Pens

The Purple Surgical Skin Marker Pens allow for precise marking on skin prior to and during surgery.

Features and benefits include:

- Sterile, single-use convenience
- Gentian violet colour ink is quick drying and does not run
- Enclosed flexible ruler ensures accurate measurement, even over contoured surfaces
- Large ink reservoir
- Full range of surgical applications due to the availability of standard, fine and dual tip options



Product	Code	Description	Pack Qty
Skin Marker Pens	PS3150	Skin Marker Standard Tip with Flexible Ruler	25 / Box
Skin Marker Pens	PS3150L	Skin Marker Standard Tip with Flexible Ruler and Label	25 / Box
Skin Marker Pens	PS3151	Skin Marker Fine Tip with Flexible Ruler	25 / Box
Skin Marker Pens	PS3151L	Skin Marker Fine Tip with Flexible Ruler and Label	25 / Box
Skin Marker Pens	PS3158	Skin Marker Dual Tip with Flexible Ruler	25 / Box
Skin Marker Pens	PS3158L	Skin Marker Dual Tip with Flexible Ruler and Label	25 / Box

Swab Counting Bags

The Purple Surgical Swab Counting Bags allow for the safe, effective, convenient counting, handling and disposal of contaminated swabs.

Features and benefits include:

- Safe, hygienic system minimises the risk of pathogen inhalation
- Clear construction allows visualisation for rapid, accurate counting
- Large or small pocket options
- Enables simple swab weight measurement for blood loss analysis
- Courtesy opening tab minimises handler's exposure to contaminants
- Available with standard or extended width swab racks




Product	Code	Description	Pack Qty
Swab Counting Bags	PS3050	Swab Bags, Non-Sterile	250 / Box
Swab Counting Bags	PS3055	Swab Rack, Reusable	1 / Each
Swab Counting Bags	PS3056	Swab Rack Extension Bar, Reusable	1 / Each

Abrasive Tip Cleaner

The Purple Surgical Abrasive Tip Cleaner allows for the effective removal of charred debris from electrodes.

Features and benefits include:

- Large 5cm x 5cm effective area
- Extended release liner enables convenient removal with gloved hands
- Radio-opaque marker ensures patient safety through rapid detection
- Coated abrasive surface ensures no flaking of contaminated abrasive material during use




Product	Code	Description	Pack Qty
Abrasive Tip Cleaner	P55200	5cm x 5cm	25 / Box
Abrasive Tip Cleaner	P55204	5cm x 5cm	100 / Box

Electrosurgical Pencils

The Purple Surgical Electrosurgical Pencils are available as hand operated, or foot operated, with a full range of accessories.

Features and benefits include:

- Sterile, single-use convenience
- Comfortable, ergonomic design
- Tactile cut and coag button action
- Universal ESU fitting
- Extensive range of electrode accessories



Product	Code	Description	Pack Qty
Electrosurgical Pencil	P55555	Handswitch Operated Pencil	25 / Box
Electrosurgical Pencil	P55550	Footswitch Operated Pencil	25 / Box
Electrosurgical Pencil Accessories	P55003	Ball Electrode	50 / Box
Electrosurgical Pencil Accessories	P55004	Extended Ball Electrode	50 / Box
Electrosurgical Pencil Accessories	P55005	Flat Blade Electrode	50 / Box
Electrosurgical Pencil Accessories	P55006	Extended Flat Blade Electrode	50 / Box
Electrosurgical Pencil Accessories	P55007	Needle Electrode	50 / Box
Electrosurgical Pencil Accessories	P55009	Teflon Coated Flat Blade Electrode	50 / Box

Light Handle Covers

The Purple Surgical Light Handle Covers are single-use, sterile covers that fit easily and securely over our universal operating theatre light handle, which is supplied on a 'free loan' basis.

Features and benefits include:

- Sterile, single-use convenience
- Gives control of light direction to surgeon and/or scrub personnel
- Securely retained by the unique 'rigid' flange and ribs
- Full range of handles ensure compatibility with all theatre lights
- Ability to conveniently 're-cover' if contaminated
- Available in single, double or triple packs for maximum convenience and minimal waste



Product	Code	Description	Pack Qty
Light Handle Covers	P54611	Single Pack	100 / Box
Light Handle Covers	P54612	Dual Pack	100 / Box
Light Handle Covers	P54613	Triple Pack	100 / Box

Vessel Loops

Purple Surgical Vessel Loops allow for the atraumatic retraction, occlusion and identification of arteries, veins, nerves and ureters.

Features and benefits include:

- Soft, strong, atraumatic silicone construction
- Oval cross section ensures even dispersion of pressure
- Radio-opaque material ensures patient safety through rapid location
- Colour coded for use with arteries, veins, nerves and ureters
- Available in two sizes – mini and maxi



Product	Code	Description	Colour	Pack Qty
Silicone Vessel Loops	P53001	Mini	White	10 / Box
Silicone Vessel Loops	P53002	Mini	Red	10 / Box
Silicone Vessel Loops	P53003	Mini	Blue	10 / Box
Silicone Vessel Loops	P53004	Mini	Yellow	10 / Box
Silicone Vessel Loops	P53201	Maxi	White	10 / Box
Silicone Vessel Loops	P53202	Maxi	Red	10 / Box
Silicone Vessel Loops	P53203	Maxi	Blue	10 / Box
Silicone Vessel Loops	P53204	Maxi	Yellow	10 / Box

Purple Arm Slings

The Purple Surgical Arm Slings are a purpose designed arm elevation sling that provides comfortable support after burns, surgery or trauma.

Features and benefits include:

- Patient comfort is assured due to the soft foam construction
- Wrist support is achieved by means of the extended lower arm section
- Velcro fastening allows simple, quick wound site inspection
- Design enables arm elevation to promote fluid drainage
- Machine washable
- Available in both adult and child sizes



Product	Code	Description	Pack Qty
Purple Arm Sling	PS9000	Adult Size	10 / Box
Purple Arm Sling	PS9001	Child Size	10 / Box

Eye Shields

The Purple Surgical Eye Shields have been designed for facial protection against blood, splashes, body fluids and chemicals.

Features and benefits include:

- Lightweight frame with flexible arms for user comfort
- Individually wrapped for single-use application
- Available pre-assembled or in two parts (lens & frame) for simple self-assembly



Product	Code	Description	Pack Qty
Eye Shield	PS9220	Pre-Assembled	20 / Box
Eye Shield	PS9210	Self-Assembly	50 / Box

Face Visor

The Purple Surgical Face Visor has been designed for full facial protection against blood, splashes, body fluids and chemicals.

Features and benefits include:

- Lightweight visor with elastic strap for secure positioning and user comfort
- Foam head band provides ventilation to prevent misting
- Available pre-assembled and individually wrapped for single-use application



Product	Code	Description	Pack Qty
Face Visor	PS9200	Pre-Assembled	25 / Box

Non-Conductive Suction Connection Tubes

Purple Surgical Suction Connection Tubes are manufactured from clear PVC, with moulded connectors. Tubes are supplied double wrapped in high density polyethylene inner bag, for simple aseptic removal, with easy peel pouch.

Features and benefits include:

- Supplied with male and female connectors
- Available in 5mm, 6mm or 7mm internal diameters
- Adhere to current REACH guidelines



Product	Code	Description	Length x Diameter	Pack Qty
Suction Connection Tubes	PSST/06/FFM5	Non-Conductive, Female / Female / Male	1.8m (6ft) x 5mm	25 / Box
Suction Connection Tubes	PSST/06/FFM6	Non-Conductive, Female / Female / Male	1.8m (6ft) x 6mm	25 / Box
Suction Connection Tubes	PSST/06/FFM7	Non-Conductive, Female / Female / Male	1.8m (6ft) x 7mm	25 / Box
Suction Connection Tubes	PSST/08/FFM5	Non-Conductive, Female / Female / Male	2.5m (8ft) x 5mm	25 / Box
Suction Connection Tubes	PSST/08/FFM6	Non-Conductive, Female / Female / Male	2.5m (8ft) x 6mm	25 / Box
Suction Connection Tubes	PSST/08/FFM7	Non-Conductive, Female / Female / Male	2.5m (8ft) x 7mm	25 / Box
Suction Connection Tubes	PSST/10/FFM5	Non-Conductive, Female / Female / Male	3m (10ft) x 5mm	25 / Box
Suction Connection Tubes	PSST/10/FFM6	Non-Conductive, Female / Female / Male	3m (10ft) x 6mm	25 / Box
Suction Connection Tubes	PSST/10/FFM7	Non-Conductive, Female / Female / Male	3m (10ft) x 7mm	25 / Box
Suction Connection Tubes	PSST/12/FFM5	Non-Conductive, Female / Female / Male	3.65m (12ft) x 5mm	25 / Box
Suction Connection Tubes	PSST/12/FFM6	Non-Conductive, Female / Female / Male	3.65m (12ft) x 6mm	25 / Box
Suction Connection Tubes	PSST/12/FFM7	Non-Conductive, Female / Female / Male	3.65m (12ft) x 7mm	25 / Box

Non-Conductive Coil Tubing

Purple Surgical Non-Conductive Coils are supplied in bulk packs with a bubble every metre. The tubing can be cut to provide a funnel or tapered connector as required. Coils are packed in a polyethylene bag within a dispenser box, with one end presented for easy extraction.

Features and benefits include:

- Available in 5mm, 6mm or 7mm internal diameters
- Supplied in 30 and 50 metre lengths
- Adhere to current REACH guidelines



Product	Code	Description	Length x Diameter	Pack Qty
Coil Tubing	PSST/30/5	Non-Sterile, Non-Conductive	30m x 5mm	1 / Each
Coil Tubing	PSST/30/6	Non-Sterile, Non-Conductive	30m x 6mm	1 / Each
Coil Tubing	PSST/30/7	Non-Sterile, Non-Conductive	30m x 7mm	1 / Each
Coil Tubing	PSST/50/5	Non-Sterile, Non-Conductive	50m x 5mm	1 / Each
Coil Tubing	PSST/50/6	Non-Sterile, Non-Conductive	50m x 6mm	1 / Each
Coil Tubing	PSST/50/7	Non-Sterile, Non-Conductive	50m x 7mm	1 / Each

Conductive Suction Connection Tubes

Purple Surgical Conductive Suction Connection Tubes are manufactured from clear PVC with a carbon based conductive line and connectors. Tubes are supplied double wrapped in a high density polyethylene inner bag, for simple aseptic removal, within an easy peel pouch.

Features and benefits include:

- Supplied with male and female connectors
- Available in 5mm, 6mm or 7mm internal diameters
- Adhere to current REACH guidelines



Product	Code	Description	Length x Diameter	Pack Qty
Suction Connection Tubes	PSCST/06/FFM5	Conductive, Female / Female / Male	1.8m (6ft) x 5mm	25 / Box
Suction Connection Tubes	PSCST/06/FFM6	Conductive, Female / Female / Male	1.8m (6ft) x 6mm	25 / Box
Suction Connection Tubes	PSCST/06/FFM7	Conductive, Female / Female / Male	1.8m (6ft) x 7mm	25 / Box
Suction Connection Tubes	PSCST/08/FFM5	Conductive, Female / Female / Male	2.5m (8ft) x 5mm	25 / Box
Suction Connection Tubes	PSCST/08/FFM6	Conductive, Female / Female / Male	2.5m (8ft) x 6mm	25 / Box
Suction Connection Tubes	PSCST/08/FFM7	Conductive, Female / Female / Male	2.5m (8ft) x 7mm	25 / Box
Suction Connection Tubes	PSCST/10/FFM5	Conductive, Female / Female / Male	3m (10ft) x 5mm	25 / Box
Suction Connection Tubes	PSCST/10/FFM6	Conductive, Female / Female / Male	3m (10ft) x 6mm	25 / Box
Suction Connection Tubes	PSCST/10/FFM7	Conductive, Female / Female / Male	3m (10ft) x 7mm	25 / Box
Suction Connection Tubes	PSCST/12/FFM5	Conductive, Female / Female / Male	3.65m (12ft) x 5mm	25 / Box
Suction Connection Tubes	PSCST/12/FFM6	Conductive, Female / Female / Male	3.65m (12ft) x 6mm	25 / Box
Suction Connection Tubes	PSCST/12/FFM7	Conductive, Female / Female / Male	3.65m (12ft) x 7mm	25 / Box

Conductive Coil Tubing

Purple Surgical Conductive Coils are supplied in bulk packs with a bubble every metre. The tubing can be cut to provide a funnel or tapered connector as required. Coils are packed in a polyethylene bag within a dispenser box, with one end presented for easy extraction.

Features and benefits include:

- Available in 5mm, 6mm or 7mm internal diameters
- Supplied in 30 and 50 metre lengths
- Adhere to current REACH guidelines



Product	Code	Description	Length x Diameter	Pack Qty
Coil Tubing	PSCST/30/5	Non-Sterile, Conductive	30m x 5mm	1 / Each
Coil Tubing	PSCST/30/6	Non-Sterile, Conductive	30m x 6mm	1 / Each
Coil Tubing	PSCST/30/7	Non-Sterile, Conductive	30m x 7mm	1 / Each
Coil Tubing	PSCST/50/5	Non-Sterile, Conductive	50m x 5mm	1 / Each
Coil Tubing	PSCST/50/6	Non-Sterile, Conductive	50m x 6mm	1 / Each
Coil Tubing	PSCST/50/7	Non-Sterile, Conductive	50m x 7mm	1 / Each

Oxygen Connecting Tubes

The Purple Surgical Oxygen Connecting Tubes are manufactured from soft green tinted PVC with a star shaped, kink restraint internal lumen and flexible moulded female connectors. This tubing is specifically designed for easy connection to oxygen therapy devices.

Features and benefits include:

- Individually wrapped
- Alternatively available with formed half bubble ends
- Available in 1.8 and 3 metre lengths
- Adhere to current REACH guidelines



Product	Code	Description	Length x Diameter	Pack Qty
Oxygen Connecting Tubes	PSOT/06/FF3	Non-Sterile, Connecting Tube	1.8m (6ft) x 3mm	50 / Box
Oxygen Connecting Tubes	PSOT/10/FF3	Non-Sterile, Connecting Tube	3m (10ft) x 3mm	50 / Box

Oxygen Multi-Fit Connecting Tubes

The Purple Surgical Multi-Fit Oxygen Connecting Tubing is manufactured from soft green tinted PVC. The tube can be cut to provide a perfect fit to oxygen therapy devices. Multi-Fit tubing is supplied in bulk coils of 30 metre or 100 metre lengths with a bubble every metre. Coils are wrapped in a polyethylene bag within a dispenser box with the tubing presented for easy extraction.

Features and benefits include:

- Available in 30 and 100 metre lengths
- Adhere to current REACH guidelines



Product	Code	Description	Length x Diameter	Pack Qty
Oxygen Multi-Fit Connecting Tubes	PSOT/30/3	Non-Sterile, Multi-Fit Connecting Tube	30m x 3mm	1 / Each
Oxygen Multi-Fit Connecting Tubes	PSOT/100/3	Non-Sterile, Multi-Fit Connecting Tube	100m x 3mm	1 / Each

